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Here is the illustrated user manual of this microtome; we hope it will be understandable.
In some time (as soon as possible), we will make an explanatory video that will be even more telling.

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The plate has in its center a hollow cylinder which circulates as a moving part, from top to bottom; it contains jaws intended to receive the specimens

to be cut.

Be careful! When handling the microtome, always hold it in an upright position, as over time the glass plate may peel off and break if it falls to the ground.

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
There are two moving parts:

A large knob graduated, in principle, in microns; but it is very relative! From experience, we know that the complete lathe allows a slice with a thickness of 10-12 microns. When you turn the knob from left to right, clockwise, the cylinder will go up.

A smaller knob to open or close the jaws that will receive the material. This knob indicates the position of the cylinder. In this case, the cylinder is completely reassembled and it will be impossible to still make cuts.

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When you have made all your cuts, the moving cylinder is completely up positioning and the jaws are flush with the glass plate (see the photo above).

To be able to make new cuts, it is necessary to perform two operations:

+++ Turn the graduated knob from right to left (counter clockwise) until you see a good part of the screw pitch (1). It is important to know that when you perform this operation, the cylinder does not automatically go down to a low position.

+++ With the finger, press down strongly on the jaws, in order to manually lower the cylinder; the position of the screw (2) indicates that the cylinder is well lowered.

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
On the picture (previous page), you see the low position of the cylinder when you pushed it down with your finger. The jaws are open. It is now time to prepare the

specimen for cutting. For this, we use extruded polystyrene (that is, sufficiently hard) about 1 cm thick. This is easily found in DIY stores. Using a cutter, cut strips 3 mm thick, and then about 2 cm long. Take the two prepared pieces and put the specimen in between.

Position the assembly in the cylinder and then close the jaws with the screw as shown in the photo above. One can clearly see between the two polystyrene plates, a blade of Basidiomycete that is colored in Congo red.

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For this, we will use a tool with a blade about 7 cm long, and is very rigid (1): it will allow to make cuts without ripples. An important precaution:

+++ Prepare a small container filled with water (for example, the lid of a 3 cm Ø Petri dish). We will deposit all the cuts ; in the water, they will spread, and we can then observe everything with a stereoscopic magnifying glass to choose the best.

+++ Using a pipette, put a drop of water on the polystyrene (2). Make the cut in a full movement from left to right or right to left depending on whether you are left-handed or right-handed, starting at the end of the blade. Above all, do not perform a saw movement.

+++ Your cup is now on the blade (3); put a drop of water on it with a pipette, so that it does not stick to the metal. Collect it with a very thin tong and put it in the container filled with water.

+++ Choose the best cut and then observe under a microscope, preferably in glycerinated water.